





SELEC

WASHROOM CLEANER

Concentrated Washroom Cleaner

Eco Select™ Concentrated Washroom Cleaner provides control to routine hygiene operations. The specially designed bottle enables exact dosing and optimised economy. Kills 99.99% of germs including Listeria, MRSA and E.Coli*.

The fresh fragranced formulation cuts through grease, grime and soap scum and removes limescale from sinks and taps leaving them gleaming and fresh. When diluted 1:60 it provides cleaning and disinfection of all hard surfaces including tiles, sinks, toilet bowls and urinals.

- Kills 99.99% of germs (passes BS EN 1276)
- Removes limescale & soap scum
- Cleans, disinfects & leaves a sparkling shine
- Forest fragrance
- Suitable for concentrated machine dispensers

(* Diluted 1:60).







100

Trigger

Spravs

BIODEGRADABILITY:



All surfactants used in Super Professional

products comply with the current

European regulations concerning

biodegradability & protection of the



50

1 Litre Makes

clockwise until tight. Directions for use: For Dosing Machines: Refer to manufacturers machine instructions for installation of bottle.

For trigger bottles: Add 600ml of clean water into the trigger bottle. Then add 1 x 10ml dose of product. Shake gently to mix. To ensure disinfection, leave for 5 minutes.

For general cleaning: Add 5 litres of warm water into a bucket. Then add 1x 20ml dose of product. Stir to mix.

Eco Select™ Concentrated Washroom provides control to routine hygiene operations. For use in manual chemical dosing systems, enabling exact dosing and optimised economy. Kills 99.99% of germs including Listeria, MRSA and E.Coli*.

The fresh fragranced formulation cuts through grease, grime and soap scum and removes limescale from sinks and taps leaving them gleaming and fresh. When diluted 1.60 if provides cleaning and disinfection of all hard surfaces including tiles, sinks, tollet bowls and urinals.

Cap instructions: To remove cap: Place bottle on a firm flat surface. Grip cap and turn anti-clockwise. Note: Use only as instructed. Do not use undiluted. * Diluted 1:60 To replace cap: Place bottle on a firm flat surface. Screw cap Storage: Store sealed in original container in a cool, dry place awa from sunlight.

Product not intended for amateur use. Safety Data Sheet available on request. This product has not been tested o animals.

BOTTLE - HDPE LABEL - PP

DANGER

clothing/eye
protection/face
protection,
FS WALLOWED: Brinse mouth,
Do NOT induce vomiting,
FO SISM (or harb) Take off
immediately all contaminated clothing,
Rinse skin with water or shower.
FI NE YES: Rinse cautiously with water for
several minutes. Remove contact lenses,
if present and easy to do. Continue
(rinsing).

rinsing.
Immediately call a POISON
CENTRE/doctor. Dispose of
contents/container in accor
national regulations.


Contains: Citric Acid Anhydrous, Alkyl (C12–16) Dimethylbenzyl Ammonium Chloride (ADBAC/BKC (C12–16)) 3.38g/100g. 2–(2–Butoxyethoxy) Ethanol. Alcohols, C12–14, Ethoxylated.

ORDER CODE:

2ltr - 800-299-0013



Box qty per case: 4 x 2Ltr Bottles



100 cases to a pallet 5 layers to a pallet

QUALITY ASSURANCE:

This product is manufactured in the UK by MIRIUS, A Coventry Group Company.

Produced under ISO 9001 Quality Management System & ISO 14001 **Environmental Management System. This** ensures our products and services are of the highest possible standard.

This product has not been tested on animals.



environment



Scan Me

mirius

Date: 15/10/2019













